# Sensation™ 7Fr.

---

## INSTRUCTIONS FOR USE

Read this entire booklet before opening inner pack

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

ENGLISH



EXPLANATION OF SYMBOLS ON PACKAGE LABELING

	Do Not Reuse
REF	Catalog Reference Number
SN	Product Serial Number
	Lot/Batch Code
	Date of Manufacture (YYYY-MM-DD)
	Expiration Date (YYYY-MM-DD)
	Sterile. Ethylene Oxide Sterilization Method
	See Instruction for Use
	CE Conformity Marking in Accordance with the European Council Directive 93/42/EEC Concerning Medical Devices

Sensation is a trademark of Datascope Corp.

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
INTENDED USE

The intra-aortic balloon catheter and accessories are used to provide counterpulsation therapy in the aorta, whereby balloon inflation during diastole and deflation during systole increases blood supply to the heart muscle and decreases the work of the left ventricle.

PRODUCT STORAGE AND HANDLING REQUIREMENTS

Product sterile if packaging is undamaged or unopened.  
DO NOT sterilize.  
Handle with care. Do not expose to direct sunlight.

I. Indications for use

- A. Refractory Unstable Angina.
- B. Impending Infarction.
- C. Acute Myocardial Infarction.
- D. Refractory Ventricular Failure.
- E. Complications of Acute MI (i.e. Acute MR or VSD or papillary muscle rupture).
- F. Cardiogenic Shock.
- G. Support for diagnostic, percutaneous revascularization and interventional procedures.
- H. Ischemia related intractable ventricular arrhythmias.
- I. Septic Shock.
- J. Intraoperative pulsatile flow generation.
- K. Weaning from cardiopulmonary bypass.
- L. Cardiac support for non-cardiac surgery.
- M. Prophylactic support in preparation for cardiac surgery.
- N. Post-surgical myocardial dysfunction/low cardiac output syndrome.
- O. Myocardial Contusion.
- P. Mechanical bridge to other assist devices.
- Q. Cardiac support following correction of anatomical defects.

II. Contraindications

- A. Severe aortic insufficiency.
- B. Abdominal or aortic aneurysm.
- C. Severe calcific aorta-iliac disease or peripheral vascular disease.
- D. Introduction of the IAB catheter without the use of an introducer sheath is not recommended in patients with severe obesity, scarring of the groin or other contraindications to percutaneous insertion.

III. Summary of warnings and precautions

A WARNINGS

- 1. If you continue to pump an IAB that has a leak, gaseous embolic injury of organs may result or a large blood clot may form within the balloon membrane requiring surgical removal of the IAB catheter.
- 2. Do not inflate the IAB using a syringe or any other means if a balloon membrane leak is suspected.
- 3. Perforation of a balloon membrane may indicate that the patient's vascular condition may induce abrasion or perforation in subsequent balloon membranes.
- 4. The practitioner must be aware of adverse effects associated with percutaneous sheath introduction including bleeding at the insertion site, limb ischemia, infection, vessel trauma, and thrombosis.
- 5. The practitioner must be aware of the potential for air embolism associated with open needles, sheaths, or catheter lumens in the patient's vasculature.
- 6. Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body fluid precautions in the care of all patients.
- 7. Do not cut the guidewire.
- 8. Use only the supplied introducer sheath, other commercially available sheaths are not recommended for the insertion of the Sensation 7Fr IAB catheter.

- 9. Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.
- 10. Use only the dilator that is packaged with the introducer set.
- 11. Use the supplied 0.035" (0.09cm) guide wire for insertion of the introducer only.
- 12. Do not cut the introducer sheath to alter the length.
- 13. Blood leakage may occur if the hemostasis valve is removed from or not secured properly to the introducer sheath.
- 14. Whenever possible, use fluoroscopy during IAB catheter insertion to ensure proper placement.
- 15. Do not insert the IAB catheter unless the inner lumen is supported by guidewire.
- 16. Do not use excessive force when inserting the IAB catheter. If you use excessive force when inserting the IAB catheter, arterial tearing, dissection, or balloon membrane damage may occur.
- 17. Any kinking or damage to the inner lumen may result in subsequent fatigue failure to the inner lumen when pumping.
- 18. If you did not use fluoroscopy during IAB catheter insertion, you MUST take an X-ray as soon as possible to ascertain that the IAB catheter is properly positioned.
- 19. If the balloon membrane has not entirely exited the introducer sheath, it will not inflate and deflate properly.
- 20. Do not insert any portion of the Universal Sheath Seal below the skin line.
- 21. NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).
- 22. Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.
- 23. If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediate removal of the IAB catheter is indicated as damage may have occurred to the balloon membrane during insertion.
- 24. If you note any restraint, unusual folding patterns of the balloon membrane or flexing of the inner lumen immediately reposition the IAB catheter. The life of a balloon membrane may be unpredictably shortened as a result of restraint, which could lead to a leak.
- 25. Do not attempt to withdraw the balloon membrane through the introducer sheath.
- 26. If you feel any undue resistance during withdrawal of the IAB catheter, discontinue withdrawal and consider removal of the IAB catheter via an arteriotomy. Difficult removal may be a result of entrapment, due to a dried blood clot having formed within the balloon membrane from a balloon membrane leak.
- 27. If limb ischemia is observed after IAB catheter removal, a vascular procedure may be indicated.

B. Precautions

- 1. Whenever possible, use fluoroscopy during guide wire and introducer sheath insertion.
- 2. Do not remove the T-handle or packaging stylet until immediately prior to insertion.
- 3. The stylet wire provides support to the IAB. Handle the IAB with care and be sure to support the t-handle so as not to kink or place undue force on the catheter.
- 4. Remove the IAB catheter from the T-handle by pulling STRAIGHT out to avoid damaging it.
- 5. Take care not to kink or place undue force on the IAB catheter.
- 6. DO NOT handle the IAB Membrane or wipe the catheter prior to insertion.
- 7. DO NOT twist the catheter during insertion.
- 8. Use the appropriately sized guidewire for the intra-aortic balloon catheter. Use only the 0.018" (0.05 cm) guide wire with the Sensation 7 Fr. IAB Catheter, for IAB catheter insertion, whether performing sheathless or sheathed insertion.
- 9. Always advance in short, continuous, one inch (2.5 cm) strokes to avoid kinking the IAB catheter.
- 10. Use care to prevent kinking of the introducer sheath during insertion.
- 11. Pinching or kinking the reinforced introducer sheath may damage it, preventing insertion of the IAB catheter.
- 12. If using a chest X-ray to identify the location of the IAB catheter, it is recommended to place the IABP in the standby mode while the X-ray is taken and then immediately resume pumping.
- 13. Do not place any sutures around the outside diameter of the introducer to avoid kinking or damaging the sheath.
- 14. If the inner lumen has been capped off DO NOT attempt to re-access the inner lumen.
- 15. Take care not to kink the optical sensor cable to avoid damage.
- 16. DO NOT touch the exposed end of the optical sensor cable, or permit it to contact other surfaces. This could damage or contaminate the sensor connection.
- 17. The inner lumen should be connected to a standard flushing apparatus. This may reduce incidence of occlusion due to blood coagulation. The use of a transducer to monitor the inner lumen is optional.
- 18. The IAB should not remain inactive (i.e., not inflating and deflating) for more than 30 minutes because of the potential for thrombus formation.
- 19. To avoid cutting the IAB catheter or the introducer, do not use scissors to remove the dressing.
- 20. Inspect the entire introducer and IAB catheter upon removal to be certain the entire device has been removed.

IV. Adverse effects

A. Balloon Membrane Perforation

- Balloon membrane perforation may be caused by:
- Contact with a sharp instrument
- Fatigue failure due to unusual (biaxial) folding of the balloon membrane during use
- Contact with calcified plaque resulting in abrasion of the surface and eventual perforation

If perforation occurs, blood may be visible in the IAB catheter. If balloon membrane perforation is suspected, as may be evidenced by: 1) IAB pump leak alarms, 2) dried blood particles or serosanguineous fluid seen in the extracorporeal tubing or catheter extender or 3) a sudden change in the diastolic augmentation pressure waveform, the following procedure must be performed immediately:

- 1. Stop pumping.
- 2. Remove the IAB catheter.
- 3. Placing the patient in the Trendelenburg position should be considered if a leak is suspected.
- 4. Consider IAB catheter replacement if the patient's condition warrants.

WARNINGS:

- If you continue to pump an IAB with a leak, gaseous embolic injury of organs may result or a large blood clot may form within the balloon membrane requiring surgical removal of the IAB catheter.
- Do not inflate the IAB using a syringe or any other means if a balloon membrane leak is suspected.
- Perforation of a balloon membrane may indicate that the patient's vascular condition may induce abrasion or perforation in subsequent balloon membranes.

NOTE: The length of time a balloon membrane can survive such contact with plaque or unusual folding is unpredictable. A leak in an IAB catheter within the bloodstream may allow gas to enter the patient's bloodstream which may result in patient injury. Large perforations are rare, therefore the small quantity of gas released is usually asymptomatic. The rate of incidence at each individual hospital may be influenced by the degree of vascular disease in that patient population, by the location of the IAB catheter in the aorta, or by using a balloon membrane size inappropriate for the specific patient.

B. Limb Ischemia

During or after IABC therapy, limb ischemia may result. It can be caused by an obstruction of flow due to:

- Thrombus formation
- Creation of an intimal layer separation or flap
- The presence of the introducer sheath or IAB catheter

After IAB catheter removal, if limb ischemia is observed, a vascular procedure may be indicated. Monitor distal limb for the development of compartment syndrome.

C. Bleeding at the Insertion Site

Bleeding at the insertion site may be caused by:

- Trauma to the artery during insertion of the IAB
- Excessive catheter movement at the insertion site
- Anticoagulation

Bleeding at the insertion site may be controlled with direct pressure at the insertion site, assuring adequate distal blood flow. If bleeding persists, surgical repair of the insertion site may be indicated.

D. Infection

Infection may occur due to interruption of normal skin integrity at the IAB catheter insertion site. Sterile technique should be used during IAB catheter insertion and during dressing changes. Assess the patient for the development of IAB catheter related infection and treat, if necessary.

E. Thrombocytopenia

Thrombocytopenia may develop due to mechanical damage to the platelets. Monitor platelet count and replace platelets, if necessary.

F. Aortic Dissection

Aortic dissection may occur during insertion of the IAB catheter. Symptoms can include back and/or abdominal pain, a decreased hematocrit, and hemodynamic instability.

G. Thrombosis

Thrombus formation may occur during counterpulsation. The symptoms associated with thrombosis formation and treatment will depend on the organ system involved.

V. Equipment required

The following is a list of compatible consoles for the Datascope Sensation 7Fr. IAB catheter at heart rates not exceeding 140 BPM. The use of the Datascope Sensation 7Fr. catheter on some consoles at higher heart rates may result in decreased augmentation.

NOTE: Optical sensor compatibility only with Datascope CS300 console. A conventional arterial pressure signal (radial artery) or external pressure signal is required to establish an arterial pressure waveform on all other pumps.

IAB Catheter	Datascope Systems	*Arrow Pumps
Sensation 7 Fr.	98, 98XT, CS100 and CS300	Arrow ACAT and AutoCat

\*Arrow/Kontron Catheter Extender, PN 0684-00-0112-01, must be ordered separately.

IAB catheter sizing

Select the most appropriate size IAB catheter for the patient from the following balloon membrane sizing table:

NOTE: This information is to be used only as a guideline. Clinical judgment and patient factors (i.e., torso length) should be considered when selecting the most appropriate size IAB catheter.

Ballon Membrane Volume (cc)	Ballon Membrane Dimensions		Patient Height	
	Length (mm)	Diameter (mm)	(ft)	(cm)
34	220	15	5'0"-5'4"	152-162
40	255	15	5'4"-6'0"	162-183

The following sterile equipment is required for insertion of the IAB catheter. Inspect all components prior to use.

The following items are supplied by Datascope and are provided with all products:

- One sterile one-way valve and 30cc syringe (included with IAB catheter)
- One sterile Insertion Kit which includes:
  - One 18 gauge angiographic needle
  - One vessel dilator
  - One introducer sheath with hemostasis valve
  - One introducer dilator
  - One three-way stopcock
  - One Datascope Catheter Extender
- Guide Wires:
  - One 0.018" (0.05cm) x 145 cm PTFE coated guide wire
  - One 0.035" (0.09 cm) x 55 cm uncoated guide wire

The following items are not supplied by Datascope:

- Local anesthetic with syringe and needle
- One sterile scalpel and blade
- One sterile 20cc syringe with saline flush
- Sterile lint free sponges
- Kelly clamp or equivalent
- One sterile 60cc syringe (optional)

## VI. Instructions

### WARNINGS:

- The practitioner must be aware of adverse effects associated with percutaneous sheath introduction including bleeding at the insertion site, limb ischemia, vessel trauma, and thrombosis.
- The practitioner must be aware of the potential for air embolism associated with open needles, sheaths, or catheter lumens in the patient's vasculature.
- Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.

**PRECAUTION:** Whenever possible, use fluoroscopy during guide wire and introducer sheath insertion.

### A. INSTRUCTIONS FOR SHEATHLESS INSERTION

**NOTE:** Insertion technique is for SHEATHLESS INSERTION of the IAB unless otherwise indicated.

**NOTE:** For SHEATHED INSERTION proceed to Section B. INSTRUCTIONS FOR SHEATHED INSERTION.

**PRECAUTION:** Do not remove the T-handle or packaging stylet until immediately prior to insertion.

1. Make the customary preparations for percutaneous catheterization and administer appropriate local anesthesia.
2. Insert the angiographic needle into the common femoral artery at a 45 degree or less angle. (See Figure 1)


Figure 1

3. Insert the J-tip end of the 0.018" (0.05cm) guide wire through the angiographic needle and advance into the thoracic aorta.

**WARNING:** Do not cut the guide wire.

**WARNING:** Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.

**PRECAUTION:** Use the appropriately sized guidewire for the insertion of the Intra-aortic balloon. Use only the 0.018" (0.05cm) guide wire with the Sensation 7 Fr. IAB Catheter, for IAB catheter insertion, whether performing sheathless or sheathed insertion.

4. Keeping the guide wire in place, remove and discard the needle.
5. Wipe the blood from the guide wire with a wet, lint-free sponge.
6. Make a small incision at the exit of the guide wire to facilitate inserting the vessel dilator through the skin.
7. Place the tapered end of the vessel dilator over the exposed guide wire and dilate the artery by pushing the vessel dilator into the arterial lumen.
8. Keeping the guide wire in place, remove and discard the vessel dilator. Apply pressure at the wound site to control bleeding.
9. Wipe the blood from the guide wire with a wet, lint-free sponge.
10. Spread the tissue at the incision with a Kelly clamp or tissue dilator.

### WARNINGS:

- Whenever possible, use fluoroscopy during IAB catheter insertion to ensure proper placement.
  - Do not insert the IAB catheter unless the inner lumen is supported by a guide wire.
11. Remove the IAB catheter tray from the sterile packaging.

**PRECAUTION:** Do not remove the T-handle or packaging stylet until immediately prior to insertion

12. Firmly attach the 30cc syringe with one-way valve to the male luer fitting of the extracorporeal tubing. (See Figure 2)


Figure 2

13. With the 30cc syringe, slowly aspirate a full 30cc (See Figure 3). Remove the syringe while leaving the one-way valve securely in place to ensure vacuum is maintained.


Figure 3

14. Carefully remove the optical sensor cable, extracorporeal tubing, Y-fitting and IAB catheter WITH T-handle from the tray – (DO NOT disconnect one-way valve when removing the extracorporeal tubing from the tray.) (See figure 4)


Figure 4

### PRECAUTIONS:

- Do not remove the T-handle or packaging stylet until immediately prior to insertion.
- Take care not to kink or place undue force on the IAB catheter.

**NOTE:** Maintain vacuum on the IAB catheter throughout insertion. Do not remove the one-way valve.

15. It is recommended that anticoagulant therapy be administered in accordance with standard hospital practice.
16. If Fluoroscopy is not used, measure the distance from the angle of Louis or between the second and third intercostal space extending down to the umbilicus and then obliquely over to the femoral insertion site. Slide the sheath seal protector up the catheter to this measured distance
17. Remove the stylet wire from the inner lumen (Figure 5). Do not attempt stylet reinsertion.

**PRECAUTION:** The stylet wire provides support to the IAB. Handle the IAB with care and be sure to support the t-handle so as not to kink or place undue force on the catheter

18. Manually flush inner lumen with 3-5cc flush solution.


Figure 5

19. Withdraw the balloon membrane from the protective T-handle by PULLING THE CATHETER STRAIGHT OUT OF THE T-HANDLE (See Figure 6)

**PRECAUTION:** Remove the IAB catheter from the T-handle by pulling STRAIGHT OUT to avoid damaging it.


Figure 6

**PRECAUTION:** DO NOT handle the IAB Membrane or wipe the catheter prior to insertion.

20. Insert the 0.018" (0.05cm) guide wire through the inner lumen (See Figure 7.) Advance the IAB catheter over the guide wire until the guide wire exits the female luer hub. Always ensure the Operator has complete control of the guidewire.


Figure 7

**WARNING:** Do not use excessive force when inserting the IAB catheter. If you use excessive force when inserting the IAB catheter, arterial tearing, dissection, or balloon membrane damage may occur.

**NOTE:** During insertion of the IAB catheter, arterial blood under pressure may run down the length of the folds in the balloon membrane and drip or be expelled under arterial pressure from the balloon membrane/catheter junction. THIS "CHANNELING" IS NOT A LEAK. As the IAB catheter is advanced the bleeding will diminish.

21. While controlling the proximal end of the guide wire, advance the IAB catheter over the guide wire into the artery. Always grasp the IAB catheter no more than one inch (2.5 cm) from the insertion site or sheath hub and advance in short continuous strokes to avoid kinking the IAB catheter while maintaining complete control of the guidewire.

### PRECAUTIONS:

- DO NOT twist the catheter during insertion.
- Always advance in short, continuous, one inch (2.5 cm) strokes to avoid kinking the IAB catheter

### WARNINGS:

- Any kinking or damage to the inner lumen may result in subsequent fatigue failure to the inner lumen when pumping.
  - Do not insert the IAB catheter unless the inner lumen is supported by a guidewire.
22. Advance the IAB catheter to the proper position in the descending thoracic aorta, with the IAB catheter tip just distal (approximately 2cm) to the left subclavian artery. (See Figure 8)


Figure 8

**WARNING:** If you do not use fluoroscopy during IAB catheter insertion, you MUST take an X-ray as soon as possible to ascertain that the IAB catheter is properly positioned.

**PRECAUTION:** If using a chest X-ray to identify the location of the IAB catheter, it is recommended to place the IABP in the standby mode while the X-ray is taken and then immediately resume pumping.

23. Maintain sterility of the exposed catheter until the proper IAB catheter position is verified.
24. After proper placement of the IAB Catheter, push the Universal Sheath Seal as close to the insertion site as possible. (See Figure 9).


Figure 9

**WARNING:** Do not insert any portion of the Universal Sheath Seal below the skin line.

### B. INSTRUCTIONS FOR SHEATHED INSERTION

1. Make the customary preparations for percutaneous catheterization and administer appropriate local anesthesia.
2. Insert the angiographic needle into the common femoral artery at a 45 degree or less angle. (See Figure 10)


Figure 10

3. Insert the J-tip end of the 0.035" (0.09 cm) guide wire through the angiographic needle and advance into the femoral artery.

**WARNING:**

- Do not cut the guide wire.
  - Do not withdraw the guide wire against the needle bevel to avoid possibly severing or damaging the guide wire.
  - Use only the dilator that is packaged with the introducer set.
  - Use the supplied 0.035" (0.09 cm) guide wire for insertion of the introducer ONLY.
  - Do not cut the introducer to alter the length.
  - Use only the supplied introducer sheath. Other commercially available sheaths are not recommended for the insertion of the Sensation 7Fr. IAB catheter.
- Keeping the guide wire in place, remove and discard the needle.
  - Wipe the blood from the guide wire with a wet, lint-free sponge.
  - Make a small incision at the exit of the guide wire to facilitate inserting the sheath introducer through the skin.
  - Insert the introducer dilator into the introducer sheath.
  - When using the sheath, insert the introducer dilator into the sheath hub and twist lock in place.
  - Place the tapered end of the introducer dilator over the exposed guide wire and advance the introducer with a rotary motion into the arterial lumen. (See Figure 11)


Figure 11

- Remove the dilator leaving the sheath in place.
- Remove the 0.035" (0.09cm) guide wire and replace with the 0.018" (0.05cm) guide wire and continue.

**WARNINGS:**

- Blood leakage may occur if the hemostasis valve is removed from or not secured properly to the introducer sheath
  - Whenever possible, use fluoroscopy during IAB catheter insertion to ensure proper placement.
  - Do not insert the IAB catheter unless the inner lumen is supported by a guide wire.
- Remove the IAB catheter tray from the sterile packaging.

**PRECAUTIONS:**

- Do not remove the T-handle or packaging styllet until immediately prior to insertion.
  - Use care to prevent kinking of the introducer sheath during insertion.
  - Pinching or kinking the reinforced introducer sheath may damage it, preventing insertion of the IAB catheter.
- Firmly attach the 30cc syringe with one-way valve to the male luer fitting of the extracorporeal tubing. (See Figure 12)


Figure 12

- With the 30cc syringe, slowly aspirate a full 30cc (See Figure 13). Remove the syringe while leaving the one-way valve securely in place to ensure vacuum is maintained.


Figure 13

- Carefully remove the optical sensor cable, extracorporeal tubing, Y-fitting and IAB catheter WITH T-handle from the tray – (DO NOT disconnect one-way valve when removing the extracorporeal tubing from the tray.) (See Figure 14)


Figure 14

**PRECAUTIONS:**

- Do not remove the T-handle or packaging styllet until immediately prior to insertion.
- Take care not to kink or place undue force on the IAB catheter.

NOTE: Maintain vacuum on the IAB catheter throughout insertion. Do not remove the one-way valve.

- It is recommended that anticoagulant therapy be administered in accordance with standard hospital practice.
- If fluoroscopy is not used, measure the distance from the angle of Louis or between the second and third intercostal space extending down to the umbilicus and then obliquely over to the femoral insertion site. Slide the sheath seal protector up the catheter to this measured distance
- Remove the styllet wire from the inner lumen (Figure 15). Do not attempt styllet reinsertion.

**PRECAUTION:** The styllet wire provides support to the IAB. Handle the IAB with care and be sure to support the t-handle so as not to kink or place undue force on the catheter.

- Manually flush inner lumen with 3-5cc flush solution.


Figure 15

**PRECAUTION:** DO NOT handle the IAB Membrane or wipe the catheter prior to insertion

- Withdraw the balloon membrane from the protective T-handle by PULLING THE CATHETER STRAIGHT OUT OF THE T-HANDLE (See Figure 16)

**PRECAUTION:** Remove the IAB catheter from the T-handle by pulling STRAIGHT OUT to avoid damaging it.


Figure 16

**PRECAUTION:** DO NOT handle the IAB Membrane or wipe the catheter prior to insertion.

- Insert the 0.018" (0.05cm) guide wire through the inner lumen (See Figure 17.) Advance the IAB catheter over the guide wire until the guide wire exits the female luer hub. Always ensure the Operator has complete control of the guidewire.


Figure 17

**WARNING:** Do not use excessive force when inserting the IAB catheter. If you use excessive force when inserting the IAB catheter, arterial tearing, dissection, or balloon membrane damage may occur.

NOTE: During insertion of the IAB catheter, arterial blood under pressure may run down the length of the folds in the balloon membrane and drip or be expelled under arterial pressure from the balloon membrane/catheter junction. THIS "CHANNELING" IS NOT A LEAK. As the IAB catheter is advanced the bleeding will diminish.

- While controlling the proximal end of the guide wire, advance the IAB catheter over the guide wire into the artery. Always grasp the IAB catheter no more than one inch (2.5 cm) from the insertion site or sheath hub and advance in short continuous strokes to avoid kinking the IAB catheter while maintaining complete control of the guidewire.

**PRECAUTIONS:**

- DO NOT twist the catheter during insertion.
- Always advance in short, continuous, one inch (2.5 cm) strokes to avoid kinking the IAB catheter

**WARNINGS:**

- Any kinking or damage to the inner lumen may result in subsequent fatigue failure to the inner lumen when pumping.
  - Do not insert the IAB catheter unless the inner lumen is supported by a guidewire.
- Advance the IAB catheter to the proper position in the descending thoracic aorta, with the IAB catheter tip just distal (approximately 2cm) to the left subclavian artery. (See Figure 18)


Figure 18

**WARNING:** If you do not use fluoroscopy during IAB catheter Insertion, you MUST take an

X-ray as soon as possible to ascertain that the IAB catheter is properly positioned

**PRECAUTION:** If using a chest X-ray to identify the location of the IAB catheter, it is recommended to place the IABP in the standby mode while the X-ray is taken and then immediately resume pumping.

- Maintain sterility of the exposed catheter until the proper IAB catheter position is verified.

**INSTRUCTIONS FOR USE OF SHEATH SEAL WITH INTRODUCER SHEATH**

**PRECAUTION:** Do not place any sutures or ligatures around the outside diameter of the Introducer to avoid kinking or damaging the introducer sheath.

Push the Universal Sheath Seal into the hub of the introducer sheath. (See Figure 19)


Figure 19

If it becomes necessary to reposition the IAB catheter, hold the sheath seal in one hand, grasp the catheter through the STAT-GARD sleeve with the other hand, and reposition under aseptic conditions. Do not attempt to reposition the IAB catheter by moving the sheath.

NOTE: If blood is seen passing the sheath seal following insertion through a sheath, disengage sheath seal from hemostasis valve.

**WHEN USING INTRODUCER SHEATH**

**WARNING:** If the balloon membrane has not entirely exited the introducer sheath, it will not inflate and deflate properly.

NOTE: Either during advancing or once the catheter is correctly positioned; confirm that the IAB catheter membrane has fully exited the sheath. (See Figure 20)


Figure 20

- The first single band from the IAB catheter tip indicates that the entire balloon membrane has exited from the 6" (15.2cm) introducer sheath/hemostasis valve and may now be inflated.
- Single bands follow at 3/4" (2cm) increments.

**C. INITIATING IAB PUMPING (IABP)**

NOTE: Do not raise the head of the bed greater than 45°.

- After positioning the IAB catheter, remove the guide wire.
- Once the catheter is in place, aspirate and discard 3cc of blood from the inner lumen and then immediately perform a manual flush using a syringe filled with 3cc to 5cc of flush solution. This will minimize the chances of stagnant blood clotting in the inner lumen.
- Using current hospital protocol, connect a standard arterial pressure flush apparatus to the hub of the inner lumen. A continuous 3cc/hour flow through the inner lumen is recommended

**WARNING:** NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).

**PRECAUTION:** If the inner lumen has been capped off DO NOT attempt to re-access the lumen.

- Release the vacuum from the IAB catheter by removing the one-way valve from the IAB catheter male luer fitting. (See Figure 21)


Figure 21

- Connect the IAB catheter to the pump.

**NOTES:**

- Ensure that all connections are leak-free.
- All catheter extenders are sterile and should only be used one time.
- Use one catheter extender when connecting the IAB catheter to the IAB pump

**When connecting the Sensation IAB to a Datascope CS300 IABP:**

- Connect the IAB catheter's male luer fitting to the female luer fitting of the catheter extender. Connect the male luer fitting of the catheter extender to the safety chamber/disk.
- Completely remove the cable wrap securing the optical sensor IAB cable.
- Carefully uncoil the optical sensor cable by loosely holding the coiled cable in one hand. With the other hand, grasp the optical sensor connector and slowly pull to uncoil the cable.

**PRECAUTION:** Take care not to kink the optical sensor cable to avoid damage.

6. Using the following steps, attach the optical sensor IAB cable connector to the IABP:  
NOTE: Always grasp the optical sensor cable connector to insert OR remove the connector from the pump as shown in Figure 22.

**PRECAUTION: Do NOT touch the exposed end of the Optical Sensor cable, or permit it to contact other surfaces. This could damage or contaminate the sensor connection.**

- a. With one hand, grasp the IAB's optical sensor cable connector as shown in Figure 22. Note that the red mark on the connector should be up.
- b. With the other hand, remove and discard the protective cover from the IAB optical sensor cable connector.


Figure 22

- c. With free hand, slide open the protective shutter located on the pump's rear panel by sliding it to the left to expose the IABP Sensor Input. Keep the shutter door open until connection is complete (Figure 23).
- d. Insert the optical sensor connector into the IABP's Sensor Input receptacle. The connector should be inserted until it "clicks". DO NOT touch the exposed end of the Optical Sensor Cable. (Figure 23)
- e. For securement connect the Sensor IAB cable to the tubing clips located on the catheter extender tubing.


Figure 23

NOTE: After initiating pumping on a Datascopie CS300 pump, if you observe the following help screen "Sensor IAB Failure" institute the following corrective actions:

- a. Confirm connection to pump
- b. Provide an alternate pressure source, either by:
  - i. Providing an alternate pressure source by connecting a conventional arterial pressure signal (radial artery) or external pressure signal to the pump.
  - ii. Replacing the IAB

NOTE: After initiating pumping on a Datascopie pump, if you observe the following help screen "Sensor IAB module failure" institute the following corrective actions:

- a. Provide an alternate pressure source, either by
  - i. Replacing the IABP for another sensor compatible pump.
  - ii. Providing an alternate pressure source by connecting a conventional arterial pressure signal (radial artery) or external pressure signal to the pump.
  - iii. Contact Datascopie Service for module repair.

#### When connecting the Sensation IAB using a conventional arterial pressure transducer to a Datascopie IABP:

- a. Connect the IAB catheter's male luer fitting to the female luer fitting of the catheter extender. Connect the male luer fitting of the catheter extender to the safety chamber/disk.
- b. Connect the arterial pressure cable to the pressure input port located on the rear of the pump. (See Figure 24)

NOTE: A conventional arterial pressure signal (radial artery) or external pressure signal is required to establish an arterial pressure waveform. The arterial pressure signal should derive from a site above the tip of the IAB catheter.

NOTE: A femoral signal is not recommended.


Figure 24

#### When connecting the Sensation IAB to an Arrow IAB pump:

- a. Connect the IAB catheter's male luer fitting to the female luer fitting of the Arrow catheter extender, and then connect the catheter extender to the appropriate system. Adjust the volume setting on the Arrow pump, according to the operating instructions, to match the IAB catheter volume.
- b. A conventional arterial pressure signal (radial artery) or external pressure signal is required to establish an arterial pressure waveform.

NOTE: Prior to using the Arrow catheter extender, remove the protective tubing from the pins.

7. Follow the IAB pump operating instructions to initiate pumping. If augmentation is not within the desired range, see Appendix A, Factors Affecting Augmentation.
8. If after a few cycles of counterpulsation, it appears that the balloon membrane is not fully open, perform the following procedure:

**WARNING: Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.**

- a. Detach the catheter extender from the IAB catheter's male luer fitting.
- b. Attach the supplied 3-way stopcock and a 60cc syringe to the IAB catheter's male luer fitting.

**WARNING: NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).**

- c. Aspirate to assure blood is not returned through the extracorporeal tubing.

**WARNING: If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediately remove the IAB catheter as damage may have occurred to the balloon membrane during insertion.**

- d. Inflate the IAB with air or helium, as follows, and IMMEDIATELY aspirate:

IAB Catheter Volume	Inflate Volume
34cc, 40cc	60cc

Remove the 3-way stopcock and syringe and reattach the IAB catheter's male luer fitting to the catheter extender and resume pumping.

9. Ensure that the balloon membrane is inflating and deflating in an unrestrained manner and not restrained as a result of lodging beneath a plaque, within a subintimal space, within the subclavian artery, aortic arch, abdominal aorta, or if the balloon membrane volume is too large for the specific patient's aorta.

**WARNING: If you note any restraint, unusual folding patterns of the balloon membrane or flexing of the inner lumen, immediately reposition the IAB catheter. The life of a balloon membrane may be unpredictably shortened as a result of restraint, which could lead to a balloon failure.**

10. If you observe unusual bleeding or subcutaneous hematoma at the insertion site, treat appropriately.
11. Evaluate peripheral pulses. If the distal pulse is not satisfactory or signs of limb ischemia are present, exercise discretion concerning the continuation of IAB pumping.
12. Restrain movement of the IAB catheter by securing the suture pads and the Y-fitting to the skin using a StatLock Securement Device (See Figure 25) or sutures.


Figure 25

13. Apply a dressing to the insertion site using sterile technique according to hospital policy.
14. Secure sensor cable by attaching to catheter extender using clips on extender.

#### D. MANUALLY INFLATING AND DEFLATING THE IAB CATHETER

**PRECAUTION: The IAB catheter should not remain inactive (i.e., not inflating and deflating) for more than 30 minutes because of the potential for thrombus formation.**

To keep the IAB catheter active in the event of pump failure, manually inflate and deflate the IAB as follows:

**WARNING: Do not perform manual inflation of the IAB with the catheter extender tubing attached to the IAB catheter.**

1. Detach the catheter extender from the IAB catheter's male luer fitting.
2. Attach a 3-way stopcock and a 60cc syringe to the IAB catheter's male luer fitting.

**WARNING: NEVER INJECT AIR INTO THE INNER LUMEN (female luer hub).**

3. Aspirate to assure blood is not returned through the extracorporeal tubing.

**WARNING: If you aspirate blood from the male luer fitting of the extracorporeal tubing, immediate removal of the IAB catheter is indicated as damage may have occurred to the balloon membrane during insertion.**

4. Inflate the IAB with 40cc air or helium and immediately aspirate. Repeat every 5 minutes while the IAB is inactive.
5. Remove the 3-way stopcock and syringe and reattach the IAB catheter's male luer fitting to the catheter extender and resume pumping.

#### E. REMOVING THE IAB CATHETER

1. Consider tapering or discontinuing anticoagulation therapy prior to removal.
2. Stop IAB pumping.
3. Disconnect the IAB catheter from the IAB pump permitting the IAB catheter to vent to atmosphere. Patient blood pressure will collapse the balloon membrane for withdrawal.

4. Remove all securement devices and/or sutures and dressings.

**PRECAUTION: To avoid cutting the IAB catheter or the introducer, do not use scissors to remove the dressing.**

5. Remove the IAB Catheter.

1. If an introducer sheath is used: Loosen the sheath seal from the hub and withdraw the IAB catheter through the introducer sheath until resistance is met.
2. Remove the IAB catheter and the introducer sheath as a unit.

**WARNING: Do not attempt to withdraw the balloon membrane through the introducer sheath.**  
**WARNING: If limb ischemia is observed after IAB catheter removal, a vascular procedure may be indicated.**

6. Apply digital pressure below the puncture site during IAB catheter removal. Allow free proximal bleeding for a few seconds, then apply pressure above the puncture site and allow a few seconds of back bleeding. Establish hemostasis to the puncture.

7. Carefully examine the limb distal to the insertion site for adequate perfusion.

**WARNING: If limb ischemia is observed after IAB catheter removal, a vascular procedure may be indicated.**

**PRECAUTION: Inspect the entire Introducer sheath and IAB catheter upon removal to be certain the entire device has been removed.**

NOTE: If you must initiate IAB pumping after the removal of the IAB catheter, you can perform a percutaneous insertion on the contralateral femoral artery. Do not use the same insertion site.

#### Appendix A: Factors affecting augmentation

If after pumping has commenced, augmentation is not within the desired range, one of the following may be indicated:

1. **Balloon membrane has not fully exited from the introducer sheath.** Pull the introducer sheath back until the balloon membrane fully exits the introducer sheath.
2. **Balloon membrane has not fully opened.** See VI. Instructions, Section C, Initiating IAB Pumping (IABP), Item 7.
3. **IAB augmentation/volume control on the IAB pump is set too low.** Adjust the IAB augmentation/volume control on the IAB pump.
4. **IAB catheter is positioned in the aortic arch, subclavian artery, or otherwise malpositioned in the aorta.** Observe the IAB catheter under fluoroscopy. If malpositioned, remove any securement device or suture that may have been applied to the Y-site and reposition the IAB catheter. Re-secure catheter.
5. **IAB catheter is positioned in a false lumen:** Use fluoroscopy to check for the appropriate intra-luminal position of the IAB catheter tip. If the IAB catheter is found to lie within a false lumen, remove the IAB catheter from the patient. Consider insertion of a new IAB catheter into the contralateral femoral artery:

In addition to the above, physiological conditions can contribute to poor augmentation. Among them are:

- Patient's mean arterial blood pressure is low.
- Patient's systemic vascular resistance is low.
- Patient's heart rate is rapid enough to compromise ventricular filling and ejection.

#### Appendix B: Limited warranty

Datascopie Corp. warrants that all its IAB catheters are free from defects in workmanship and materials for a period of three years from the date of purchase or until the expiration date is reached or until the date of use, whichever comes first. Datascopie Corp. shall not be liable for any incidental, special or consequential loss, damage, or expense directly arising from the use of this product. Liability under this warranty and the buyer's exclusive remedy is limited to replacement of the product which, under normal use and service, shall have been found by the Company to be defective in materials or workmanship. It shall be the buyer's obligation to return any such product to the Company for examination for replacement liability.

No agent, employee, or representative of Datascopie Corp. has any authority to bind Datascopie Corp. to any affirmation, representation, or warranty concerning its products, and any affirmation, representation, or warranty made by any agent, employee or representative shall not be enforceable by the buyer.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS, AND OF ANY OTHER OBLIGATION ON THE PART OF THE SELLER:

Damage to any product or parts through misuse, neglect, accident, or by affixing any nonstandard accessory attachments or by any customer modification voids this warranty. Datascopie Corp. makes no warranty whatever in regard to trade accessories, such being subject to the warranty of their respective manufacturers.

A condition of this warranty is that this equipment or any accessories which are claimed to be defective be returned, when authorized by Datascopie, freight prepaid, to Datascopie Corp., 15 Law Drive, CN 40011, Fairfield, New Jersey 07004-0011. Datascopie Corp. shall not have any responsibility in the event of loss or damage in transit.

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Datascope Medical Co. Ltd.
Datascope GmbH
Datascope Italia S.r.l
CAUTION: Federal (U.S.A.) Law restricts this device to sale,